

## Application Guide



**Step 1** – Find a clean soft surface such as your lawn, a large rubber mat, or a piece of carpet to boot up on, so if you select the wrong size you can exchange them unmarked. Old Macs should be fitted as the boots are marked, “left” and “right,” with the securing buckle placed on the outside of the hoof.



**Step 2** – Fully open the cover ensuring that no hook and loop areas are touching and that the cover is pulled well forward.



**Step 3** – Facing the rear of your horse, pick up the hoof and place it toe first into the boot. Pull the soft leather rear of the boot up over the bulbs of the hoof. Then, using the palm of your hand, “slap” the heel of the boot on the sole. You should feel the heels of the hoof “pop” into the boot. It should have the same feel as when you slip your own heels down into your own riding boots. If the new boot just falls onto the hoof, it is likely you have a size too big. If the heel of the hoof struggles to slip into the boot, then it is likely you have a size too small, and the heels of the hoof are sitting on the wall of the boot and not “in” the boot. Now place the hoof onto the ground and continue booting up.



**Step 4** – Wrap the main webbing strap around the pastern bone from the back of the hoof around the front and then to the back again. Feed the strap through the metal loop and then apply minimum pressure to tighten the hook and loop strap back over on itself. To confirm correct tightness of this strap, you should be able to run your finger between the strap and the pastern at the side of the pastern. It is very important that this strap is no tighter than your watch band around your wrist; otherwise it will chafe the pastern. It is also important to note here that it is the position of the cross-over on the top of the bulbs which secures the boot, not the “tightness” of it.



Once you have secured this strap, you should have a minimum of 3 inches (80mm) passed through the metal loop and folded over on itself. Any less and your boots are too small.

**Step 5** – Ensuring both side flaps are pushed firmly against the hoof wall, force the cover towards the rear of the hoof in a downward direction before crossing the ends of the cover at the rear and securing. At this point you should ensure that the foam-filled leather collar of the cover is molded firmly around the pastern bone. Also check to ensure none of the hook and loop strap is exposed at the sides. If it can be seen, the cover has been secured too high and will increase the chances of chafing at the rear of the pastern.

**Step 6** – Now secure the outside buckle strap firmly – not girth tight. Your Old Macs are now securely in place.

After you have placed your horse's hoof into the boot, you need to ensure that the entire horny hoof is sitting completely down into the boot and not “standing” on the rear of the boot. If the hoof is not sitting down comfortably onto the sole of the boot, then you have selected a size too small and need to go up a size.

### To ensure your boots are not too small:

The pastern strap, when secured, should not be too short. As previously mentioned, you should have a minimum of 3” (80mm) of the hook and loop strap folded over on itself. If you have 2” (60mm) passed back through the metal loop, you are likely to be one size too small. If you have less than 1 ½” (40mm), you are probably two sizes too small.

### To ensure your boots are not too big:

Do the “Wiggle Test.” Facing the rear of your horse, pick up the booted hoof. With the sole facing up, place a hand each side of the sole of the boot. Without using force, gently wiggle the boot clockwise and counter-clockwise. You should feel the hoof wanting to turn with the boot. If you feel the boot turning without the hoof, you have selected a size too big. You need a snug fit.

After securing the pastern strap, you should be able to place a finger between the pastern and strap at the side of the pastern bone. If it is too tight, it will chafe. Remember, it is the position of the strap which prevents the boot coming off, not the tightness.

After the cover has been secured, ensure that the hook and loop strap is not visible at the sides, and the rear ends of the cover are not secured too high. The entire cover should have no creases or folds, and should be of even height on both sides with the foam-filled leather collar snug around the pastern.


**IMPORTANT**-- It is recommended that you break your horse's new boots in slowly to reduce the chances of chafing. Old Mac's™ are no different from human leather boots or shoes which need time to soften and mold to the shape of our feet. We recommend your first few rides should be no longer than an hour, or you can turn your horse out into its paddock with the boots on for a few hours at a time and let your horse run around at its own leisure. This will allow the boots time to soften. If you know your horse has very sensitive skin and/or you are planning on riding long distances, then we recommend you use Old Mac's Gaiters.

## Gaiter Fitting Guide

Gaiters are recommended during the break in period to avoid chafing. They can be used as extra protection for sensitive horses and when riding through terrain where debris might get inside the boot.

**Step 1** – Fold the hoof plate to a 90° angle at the fold-line (see diagram), towards the inside of the Gaiter (the side without the logo).


**Steps 2-3** – Insert the Gaiter into the Old Mac's Horse Boot prior to booting the horse. Press the hoof plate down snugly onto the sole of the boot.



**Steps 4-5** – Fold the neoprene “upper” cuff down over the outside back of the boot and slip the Old Mac's Boot onto the hoof, holding the Gaiter in place with one hand. Let the horse stand down, resting the hoof flat on the ground.

**Step 6** – Unfold the Gaiter so the cuff is standing loosely upright, around the pastern.

**Step 7** – Fasten the Gaiter snugly around the pastern.





## EasyCare Inc. Warranty

EasyCare offers a 90 day repair or replacement warranty for all hoof boots (purchased as new) manufactured under the EasyCare label. Please register your product at [http://www.easycareinc.com/warranty\\_reg.aspx](http://www.easycareinc.com/warranty_reg.aspx) and keep your original dated receipt(s) as proof of purchase during the warranty period.

If a boot fails during warranty, it is generally a minor hardware issue. Please bear in mind that most components of EasyCare hoof boots are either replaceable or repairable. This includes wire laces, cables, buckles, gears, rollers, screws and straps. Most repairs can be easily accomplished with simple household tools. Depending on the circumstances, EasyCare will send you repair parts, or ask that you return the boots to us for evaluation. Videos of many repairs are available on our website at [www.easycareinc.com](http://www.easycareinc.com).

While this warranty is very broad in nature, it does not cover wear (including sole breakthrough when worn over horseshoes), misuse, abuse, factory seconds, bargain bin items, or boots purchased as used. Products returned for evaluation must be clean, dry and free of debris, dirt, sand and manure. *Failure to return clean products will result in product being shipped back at sender's expense or a cleaning fee may be assessed.* Damage inflicted by omnivorous horses, rodents or other unusual abuse is not covered.

Returns must be authorized before sending and must be returned with an RMA number. Contact EasyCare directly for all warranty issues at 1-800-447-8836 to receive an RMA number.

Specifically:

1. Returns for credit only (refund of purchase price only): Products must be in new, unused condition, free of dirt and wear marks. It is your responsibility to evaluate the fit of your new product in a manner where you may return it in new condition if you have doubts about the fit. Be certain that you know the return policy of product purchased from dealers other than direct from EasyCare (the manufacturer). EasyCare will not give credit for product purchased elsewhere. Return to point of purchase in that circumstance. Also, all brochures, inserts and/or free gifts must accompany returned items.

Online Warranty Registration:  
[easycareinc.com/warranty\\_reg.aspx](http://easycareinc.com/warranty_reg.aspx)

2. Exchanges of new product (size or color issue): EasyCare will only exchange NEW items with same product style and model for same product when requested with NEW items only. You may exchange for a different size or color exchange only, no exceptions. If you wish to have a different style of hoof boot, etc., you must return new item for refund (original point of purchase) and then purchase new choice.

3. Warranty failure (used product that is within 90 days of purchase): Contact us at 800-447-8836. EasyCare will either provide repair parts directly to you, or repair/replace after our evaluation. Items must be returned to EasyCare at the consumers cost with original receipt. EasyCare will return the items via ground service (UPS or FEDX) free of charge. Requests for faster shipping service will be at the expense of the consumer.

4. Not under warranty (over 90 days or no proof of purchase): Repair parts may be purchased from your local dealer or EasyCare, or product may be returned for evaluation and possible repair (parts, labor and freight charges will apply).

5. Call for an RMA number (under warranty or not): Returns must be authorized before sending and must be returned with an RMA number clearly written on outside of packaging. Contact us at 800-447-8836 to receive an RMA number. Returns will be refused without an RMA number on outside of package.

6. Freight charges (outside warranty period): Freight charges to and from EasyCare for repairs, exchanges and items outside the warranty period will be the purchaser's responsibility.

7. EasyCare staff is always willing to discuss any specific issues you may have regarding the performance of your EasyCare or Old Mac's products. Please contact us at 800-447-8836 or email to [admin@easycareinc.com](mailto:admin@easycareinc.com). Visit our website, [easycareinc.com](http://easycareinc.com), for tips and product information that will help with the use of any of our product line.



## Easycareinc.com 800.447.8836

Additional help is only a click away! Visit these links to learn more.

**1. Find a dealer in your area**

[easycareinc.com/Search/Dealer.aspx](http://easycareinc.com/Search/Dealer.aspx)

**2. Find a Hoof Care Professional**

[easycareinc.com/Search/Practitioner.aspx](http://easycareinc.com/Search/Practitioner.aspx)

**3. Get expert hoof boot fitting advice and the inside scoop from our blogs**

[easycareblogs.com](http://easycareblogs.com)

**4. Read the latest in hoof boot and hoof care articles**

[easycareinc.com/education/articles.aspx](http://easycareinc.com/education/articles.aspx)

**5. Check out our instructional videos**

[easycareinc.com/education/videos.aspx](http://easycareinc.com/education/videos.aspx)

**6. Subscribe to our free monthly e-newsletter**

[easycareinc.com/NewsRoom/Archives.aspx](http://easycareinc.com/NewsRoom/Archives.aspx)

**7. Download product brochures, manuals and instructions**

[easycareinc.com/education/downloads.aspx](http://easycareinc.com/education/downloads.aspx)

**8. Check out our contests and promotions**

[easycareinc.com/Cool\\_Stuff/promotions.aspx](http://easycareinc.com/Cool_Stuff/promotions.aspx)

**9. Read our press releases for breaking news**

[easycareinc.com/NewsRoom/Media.aspx](http://easycareinc.com/NewsRoom/Media.aspx)

**10. Get discount prices at our Bargain Bin**

[easycareinc.com/externallink.asp?C=XBAR](http://easycareinc.com/externallink.asp?C=XBAR)

**11. Trade in boots at the Hoof Boot Upgrade Program**

[hoofbootupgrades.com](http://hoofbootupgrades.com)

**12. Team Easyboot**

[teameasyboot.com](http://teameasyboot.com)


**13. Learn how to do your own hoof boot repairs**

[easycarererepairs.com](http://easycarererepairs.com)

We continue to strive for the best products on the market with you and your horse in mind. We hope you enjoy our product selection and website improvements. Thank you for your continued support and patronage.

**The Fine Print:**

EasyCare Inc. hereby disclaims responsibility for the misuse of any products described in this brochure and web site and/or losses incurred through the use of products by any method other than that which is outlined in this catalog or web site or EasyCare's product packaging. All products should be used in accordance with the instructions printed on the product's label. EasyCare Inc products are not intended to treat, cure or diagnose any medical condition.



## Comfort Pads

Give your horse more comfort in each step!

Comfort Pads are available in two thicknesses to allow you to decide what level of extra cushioning you give your horse. The pads are available in the original 12mm thickness, as well as the newer, thinner 6mm thickness. This thinner pad option may be used by itself or stacked in any desired combination for your custom application and is great for modifying boot fit.



“A must in EVERY boot! A perfect name for the perfect product! My navicular horse LOVED his boots with pads, they promoted healing and gave comfort no matter where we rode!

—Andy Dane

I absolutely love the comfort pads. They hold up very good and they offer my horse the frog stimulation she should have. They can also help you to get your boots to fit a little better when you just need a little extra help with the fit.

—Lisa Grigaitis

My horses love these pads. When I put on their boots with these pads in they really lick their lips when I set their foot down. After a long ride using a new set of pads, it is obvious looking at the pads how much support they have given to the horse's foot, as well as everything else up the leg. Good job EasyCare!

—Kerry Grear

Have Questions? Need Help? Visit [easycareinc.com](http://easycareinc.com)



**EasyCare**  
[easycareinc.com](http://easycareinc.com)



**G2**

Application Guide